

Rabbit Loading Control Panel Kit

Introduction

The kit contains three commonly used internal control antibodies, 30ul each.

| | Cat.# | Product | Source | Size | Application | Reactivity |
|---|--------|-----------------|--------|--------|---------------|-------------------------------------|
| 1 | AF7018 | Beta-actin Ab | Rabbit | 30ulx1 | WB,IHC,IF/ICC | Refer to the manual |
| 2 | AF7011 | Beta-tubulin Ab | Rabbit | 30ulx1 | WB,IHC,IF/ICC | Refer to the manual |
| 3 | AF7021 | GAPDH Ab | Rabbit | 30ulx1 | WB,IHC,IF/ICC | Refer to the manual |

Important Product Information

Details please refer to the product manual.

Storage

-20°C