

Goat Anti-Rabbit IgG (H+L) AP

Cat.#: S0003
Size: 100ul,200ul

Concn.: 1mg/ml
Source: Goat

Mol.Wt.:
Clonality: Polyclonal

Application: WB: 1:5000~1:10000, IHC 1:50-1:200
*The optimal dilutions should be determined by the end user.

Reactivity: Rabbit

Purification: affinity purification.

Immunogen: The antiserum was developed in goat using the Rabbit IgG as the immunogen.

Storage Condition and Buffer: AP-labeled Goat IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl.

IMPORTANT: For western blot, incubate membrane with diluted secondary Ab in 5% w/v milk , 1X TBS, 0.1% Tween@20 at 37°C, 2H.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.