

GST-Tag Antibody

Cat.#: T0007

Concn.: 1mg/ml

Mol.Wt.:

Size: 1ml,200ul,100ul,50ul

Source: Mouse

Clonality: Monoclonal

Application:

WB: 1:3000-1:10000

*The optimal dilutions should be determined by the end user.

Reactivity:

All

Purification:

Affinity-chromatography.

Specificity:


GST-Tag Mouse Monoclonal antibody detects endogenous levels of C-terminal, internal, and N-terminal GST-tagged proteins.

Immunogen:

Full-length GST protein.

Storage Condition and Buffer:

Mouse IgG1 in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.Store at -20 °C.Stable for 12 months from date of receipt.



Western blot analysis of GST-tag fusion protein, using GST-Tag Mouse Monoclonal antibody (lane 1: 1:2000; lane 2: 1:5000).

IMPORTANT: For western blot, incubate membrane with diluted primary Ab in 5% w/v milk, 1X TBS, 0.1% Tween@20 at 4°C with gentle shaking, overnight.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.