


HPLC REPORT

Sample: S-16
 Analyst: HXH
 Lot. No.: P140303-CJ566055
 Column: Gemini-NX 5 μ C18 110A, 4.6*250mm
 Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
 Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
 Gradient:

Analyzed date: 2014-03-14

	A	B
0.0min	30%	70%
25.0min	55%	45%
25.1min	100%	0%
30.0min	Stop	

Volume: 20 μ l
 Wavelength: 220nm
 Flow rate: 1.0ml/min



Rank	Time	Conc.	Area	Height
1	11.010	0.4290	27018	2515
2	11.068	0.5293	33331	4122
3	11.245	1.0792	67966	5121
4	11.632	5.3631	337748	23221
5	11.890	88.6146	5580597	646769
6	12.068	1.3953	87871	33079
7	12.150	2.2714	143043	10511
8	12.507	0.3181	20036	3055
Total		100	6297610	728393