

## ABCA1 Ab

[References\(2\)](#) [Images\(5\)](#)

Cat.#: DF8233                      Concn.: ~1mg/ml                      Mol.Wt.: 254 kDa  
Size:                                      Source: Rabbit                              Clonality: Polyclonal

Application:                      WB 1:1000-3000  
  \*The optimal dilutions should be determined by the end user.

Reactivity:                      Human,Mouse

Storage:                              Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt.

Purification:                      The antiserum was purified by peptide affinity chromatography using SulfoLink™ Coupling Resin (Thermo Fisher Scientific).

Immunogen:                      A synthesized peptide derived from human ABCA1, corresponding to a region within C-terminal amino acids.

Uniprot:                              O95477



Western blot analysis of extracts from COLO205 cells(heat-shock treatment), using ABCA1 Ab at 1/1000 dilution.  
Observed bands:250kD.

**IMPORTANT:** For western blot, incubate membrane with diluted primary Ab in 5% w/v milk , 1X TBS, 0.1% Tween@20 at 4°C with gentle shaking, overnight.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.