

## OAZ3 Ab

[Images\(3\)](#)

Cat.#: DF9726                      Concn.: ~1mg/ml                      Mol.Wt.: 27 kDa, 20 kDa  
Size:                                      Source: Rabbit                              Clonality: Polyclonal

Application:                      WB 1:1000-3000  
  \*The optimal dilutions should be determined by the end user.


Reactivity:                      Human,Mouse,Rat

Storage:                              Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt.

Purification:                      The antiserum was purified by peptide affinity chromatography using SulfoLink™ Coupling Resin (Thermo Fisher Scientific).

Immunogen:                      A synthesized peptide derived from human OAZ3, corresponding to a region within the internal amino acids.

Uniprot:                              Q9UMX2



Western blot analysis of extracts from Rat muscle, using OAZ3 Ab. The lane on the left was treated with blocking peptide.

Observed bands: 20 kDa.

**IMPORTANT:** For western blot, incubate membrane with diluted primary Ab in 5% w/v milk , 1X TBS, 0.1% Tween@20 at 4°C with gentle shaking, overnight.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.