

Goat anti-rat IgG (H+L), HRP-linked Ab

[References\(2\)](#)

Cat.#: S0009	Concn.: 1mg/ml	Mol.Wt.:
Size:	Source: Goat	Clonality: Monoclonal

Application:	WB: 1:5000~1:10000 *The optimal dilutions should be determined by the end user.
Reactivity:	Rat
Storage:	IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl. Store at -20??avoid repeated thawing and freezing.
Purification:	affinity purification.
Immunogen:	The antiserum was developed in goat using the rat IgG as the immunogen.
Conjugate:	HRP

IMPORTANT: For western blot, incubate membrane with diluted secondary Ab in 5% w/v milk , 1X TBS, 0.1% Tween@20 at 37°C, 2H.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.