

Rabbit Anti-Goat IgG (H+L) HRP

[References\(3\)](#)

Cat.#: S0010	Concn.: 1mg/ml	Mol.Wt.:
Size:	Source: Rabbit	Clonality: Polyclonal
Application:	WB: 1:5000~1:10000, IHC 1:200 *The optimal dilutions should be determined by the end user.	
Reactivity:	Goat	
Storage:	IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl. Store at -20??avoid repeated thawing and freezing.	
Purification:	affinity purification.	
Immunogen:	The antiserum was developed in Rabbit using the Goat IgG as the immunogen.	
Conjugate:	HRP	

IMPORTANT: For western blot, incubate membrane with diluted secondary Ab in 5% w/v milk , 1X TBS, 0.1% Tween@20 at 37°C, 2H.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.