

Goat Anti-Rabbit IgG (H+L) CY3-conjugated

[References\(8\)](#) [Images\(2\)](#)

Cat.#: S0011	Concn.: 1mg/ml	Mol.Wt.:
Size:	Source: Goat	Clonality: Polyclonal
Application:	WB 1:100 – 1:500, IF/ICC 1:100-1:500 *The optimal dilutions should be determined by the end user.	
Reactivity:	Rabbit	
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6.Store at -20°?avoid repeated thawing and freezing.	
Purification:	Affinity-chromatography.	
Immunogen:	The antiserum was developed in goat using the Rabbit IgG as the immunogen.	
Conjugate:	CY3	

IMPORTANT: For western blot, incubate membrane with diluted secondary Ab in 5% w/v milk , 1X TBS, 0.1% Tween@20 at 37°C, 2H.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.