

## EYFP-Tag Ab

[Images\(1\)](#)

Cat.#: BF0725                      Concn.: ~1mg/ml                      Mol.Wt.:  
Size:                                      Source: Mouse                      Clonality: Monoclonal

Application:                      WB 1:1000-1:10000  
\*The optimal dilutions should be determined by the end user.

Reactivity:                      All

Storage:                      Mouse IgG1 in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt.

Purification:                      Affinity-chromatography.



**IMPORTANT:** For western blot, incubate membrane with diluted primary Ab in 5% w/v milk, 1X TBS, 0.1% Tween@20 at 4°C with gentle shaking, overnight.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.