

## Phospho-HIPK2 (Thr273) Ab

Cat.#: AF3537  
Size:

Concn.: ~1mg/ml  
Source: Rabbit

Mol.Wt.: 140kD  
Clonality: Polyclonal

**Application:** ELISA(peptide) 1:20000-1:40000  
\*The optimal dilutions should be determined by the end user.

**Reactivity:** Human,Mouse,Rat

**Storage:** Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt.

**Purification:** The Ab is from purified rabbit serum by affinity purification via sequential chromatography on phospho-peptide and non-phospho-peptide affinity columns.

**Immunogen:** A synthesized peptide derived from human HIPK2 around the phosphorylation site of Thr273.

**Uniprot:** Q9H2X6

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.