

Phospho-HIPK2 (Tyr354) Ab

Cat.#: AF3538	Concn.: ~1mg/ml	Mol.Wt.: 140kD
Size:	Source: Rabbit	Clonality: Polyclonal
Application:	ELISA(peptide) 1:20000-1:40000 *The optimal dilutions should be determined by the end user.	
Reactivity:	Human,Mouse,Rat	
Storage:	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt.	
Purification:	The Ab is from purified rabbit serum by affinity purification via sequential chromatography on phospho-peptide and non-phospho-peptide affinity columns.	
Immunogen:	A synthesized peptide derived from human HIPK2 around the phosphorylation site of Tyr354.	
Uniprot:	Q9H2X6	

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.