

Phospho-MLK1/2 (Thr312/Thr266) Ab

Cat.#: AF3630
Size:

Concn.: ~1mg/ml
Source: Rabbit

Mol.Wt.:
Clonality: Polyclonal

Application:

ELISA(peptide) 1:20000-1:40000
*The optimal dilutions should be determined by the end user.

Reactivity:

Human,Mouse,Rat

Storage:

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt.

Purification:

The Ab is from purified rabbit serum by affinity purification via sequential chromatography on phospho-peptide and non-phospho-peptide affinity columns.

Immunogen:

A synthesized peptide derived from human MLK1/2 around the phosphorylation site of Thr312/266.

Uniprot:

P80192/Q02779

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.