

Phospho-GRK2 (Ser685) Ab

| | | |
|---------------|--|-----------------------|
| Cat.#: AF3696 | Concn.: ~1mg/ml | Mol.Wt.: |
| Size: | Source: Rabbit | Clonality: Polyclonal |
| Application: | ELISA(peptide) 1:20000-1:40000 *The optimal dilutions should be determined by the end user. | |
| Reactivity: | Human,Mouse,Rat | |
| Storage: | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt. | |
| Purification: | The Ab is from purified rabbit serum by affinity purification via sequential chromatography on phospho-peptide and non-phospho-peptide affinity columns. | |
| Immunogen: | A synthesized peptide derived from human GRK2 around the phosphorylation site of Ser685. | |
| Uniprot: | P25098 | |

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.