

Phospho-PCNA (Tyr211) Ab

Cat.#: AF3726
Size:

Concn.: ~1mg/ml
Source: Rabbit

Mol.Wt.: 36kD
Clonality: Polyclonal

| | |
|---------------|--|
| Application: | ELISA(peptide) 1:20000-1:40000 *The optimal dilutions should be determined by the end user. |
| Reactivity: | Human,Mouse,Rat |
| Storage: | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt. |
| Purification: | The Ab is from purified rabbit serum by affinity purification via sequential chromatography on phospho-peptide and non-phospho-peptide affinity columns. |
| Immunogen: | A synthesized peptide derived from human PCNA around the phosphorylation site of Tyr211. |
| Uniprot: | P12004 |

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.