

## Phospho-WEE1 (Ser642) Ab

Cat.#: AF3814	Concn.: ~1mg/ml	Mol.Wt.: 55kD(isoform), 72~100kD
Size:	Source: Rabbit	Clonality: Polyclonal
Application:	ELISA(peptide) 1:20000-1:40000 *The optimal dilutions should be determined by the end user.	
Reactivity:	Human,Mouse,Rat	
Storage:	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt.	
Purification:	The Ab is from purified rabbit serum by affinity purification via sequential chromatography on phospho-peptide and non-phospho-peptide affinity columns.	
Immunogen:	A synthesized peptide derived from human WEE1 around the phosphorylation site of Ser642.	
Uniprot:	P30291	

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.