

## Phospho-RIP3 (Thr231/Ser232) Ab

Cat.#: AF3900	Concn.: ~1mg/ml	Mol.Wt.:
Size:	Source: Rabbit	Clonality: Polyclonal
Application:	ELISA(peptide) 1:20000-1:40000 *The optimal dilutions should be determined by the end user.	
Reactivity:	Mouse,Rat	
Storage:	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt.	
Purification:	The Ab is from purified rabbit serum by affinity purification via sequential chromatography on phospho-peptide and non-phospho-peptide affinity columns.	
Immunogen:	A synthesized peptide derived from mouse RIP3 around the phosphorylation site of Thr231/Ser232.	
Uniprot:	Q9QZL0	

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.