

HIPK2 Ab

Cat.#: AF6538

Concn.: ~1mg/ml

Mol.Wt.:

Size:

Source: Rabbit

Clonality: Polyclonal

Application:

ELISA(peptide) 1:20000-1:40000

*The optimal dilutions should be determined by the end user.

Reactivity:

Human,Mouse,Rat

Storage:

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt.

Purification:

The antiserum was purified by peptide affinity chromatography using SulfoLink™ Coupling Resin (Thermo Fisher Scientific).

Immunogen:

A synthesized peptide derived from human HIPK2.

Uniprot:

Q9H2X6

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.