

## Tryptophan Hydroxylase Ab

Cat.#: AF6714  
Size:

Concn.: ~1mg/ml  
Source: Rabbit

Mol.Wt.: 55kD  
Clonality: Polyclonal

Application:

ELISA(peptide) 1:20000-1:40000

\*The optimal dilutions should be determined by the end user.

Reactivity:

Human,Mouse,Rat

Storage:

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt.

Purification:

The antiserum was purified by peptide affinity chromatography using SulfoLink™ Coupling Resin (Thermo Fisher Scientific).

Immunogen:

A synthesized peptide derived from human Tryptophan Hydroxylase.

Uniprot:

P17752

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.