

C-RAF Ab

| | | |
|---------------|--|-----------------------|
| Cat.#: AF6924 | Concn.: ~1mg/ml | Mol.Wt.: |
| Size: | Source: Rabbit | Clonality: Polyclonal |
| Application: | ELISA(peptide) 1:20000-1:40000 *The optimal dilutions should be determined by the end user. | |
| Reactivity: | Human | |
| Storage: | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt. | |
| Purification: | The antiserum was purified by peptide affinity chromatography using SulfoLink™ Coupling Resin (Thermo Fisher Scientific). | |
| Immunogen: | A synthesized peptide derived from human C-RAF. | |
| Uniprot: | P04049 | |

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.