

---

## Phospho-Phospholamban (Ser16) Ab

Cat.#: AF7458  
Size:

Concn.: ~1mg/ml  
Source: Rabbit

Mol.Wt.: 24kD  
Clonality: Polyclonal

Application: ELISA(peptide) 1:20000-1:40000

\*The optimal dilutions should be determined by the end user.

Reactivity: Human, Mouse, Rat

Storage: Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt.

Purification: The Ab is from purified rabbit serum by affinity purification via sequential chromatography on phospho-peptide and non-phospho-peptide affinity columns.

Immunogen: A synthesized peptide derived from human Phospholamban around the phosphorylation site of Ser16.

Uniprot: P26678

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.