


## hnRNP G Ab

[Images\(1\)](#)

Cat.#: DF3171                                      Concn.: ~1mg/ml                                      Mol.Wt.: 40 KD  
Size:    Source: Rabbit    Clonality: Polyclonal

Application:                                      WB 1:500-1:1000, IF/ICC 1:100-1:500  
\*The optimal dilutions should be determined by the end user.  
Reactivity:                                      Human,Mouse,Rat  
Storage:    Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt.  
Purification:                                      The antiserum was purified by peptide affinity chromatography using SulfoLink™ Coupling Resin (Thermo Fisher Scientific).  
Immunogen:                                      A synthesized peptide derived from human hnRNP G, corresponding to a region within N-terminal amino acids.  
Uniprot:    P38159



Western blot analysis of extracts from HeLa cells, using hnRNP G Ab.

**IMPORTANT:** For western blot, incubate membrane with diluted primary Ab in 5% w/v milk , 1X TBS, 0.1% Tween@20 at 4°C with gentle shaking, overnight.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.