

## FUK Ab

[Images\(1\)](#)

Cat.#: DF3347                      Concn.: ~1mg/ml                      Mol.Wt.: 110 KD  
Size:                                      Source: Rabbit                              Clonality: Polyclonal

**Application:**                      WB 1:500-1:1000, IF/ICC 1:100-1:500  
\*The optimal dilutions should be determined by the end user.


**Reactivity:**                      Human,Mouse

**Storage:**                              Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt.

**Purification:**                      The antiserum was purified by peptide affinity chromatography using SulfoLink™ Coupling Resin (Thermo Fisher Scientific).

**Immunogen:**                      A synthesized peptide derived from human FUK, corresponding to a region within N-terminal amino acids.

**Uniprot:**                              Q8N0W3



Western blot analysis of extracts from COS-7 cells, treated with insulin (0.01U/ml, 45mins), using FUK Ab.

**IMPORTANT:** For western blot, incubate membrane with diluted primary Ab in 5% w/v milk , 1X TBS, 0.1% Tween@20 at 4°C with gentle shaking, overnight.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.