

GST-Tag Ab

[References\(7\)](#) [Images\(4\)](#)

Cat.#: T0007	Concn.: ~1mg/ml	Mol.Wt.:
Size:	Source: Mouse	Clonality: Monoclonal
Application:	WB 1:3000-1:10000 *The optimal dilutions should be determined by the end user.	
Reactivity:	All	
Storage:	Mouse IgG1 in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt.	
Purification:	Affinity-chromatography.	
Immunogen:	Full-length GST protein.	

IMPORTANT: For western blot, incubate membrane with diluted primary Ab in 5% w/v milk , 1X TBS, 0.1% Tween@20 at 4°C with gentle shaking, overnight.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.