

HNF-6 Ab

[Images\(2\)](#)

Cat.#: DF2478 Concn.: ~1mg/ml Mol.Wt.: 51 kDa
Size: Source: Rabbit Clonality: Polyclonal

Application: WB 1:500-1:2000
 *The optimal dilutions should be determined by the end user.

Reactivity: Human,Mouse,Rat


Storage: Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt.

Purification: The antiserum was purified by peptide affinity chromatography using SulfoLink™ Coupling Resin (Thermo Fisher Scientific).

Immunogen: A synthesized peptide derived from human HNF-6, corresponding to a region within N-terminal amino acids.

Uniprot: Q9UBC0

Description: Transcriptional activator. Binds the consensus sequence 5'-DHWATTGAYTWW-3' on a variety of gene promoters such as those of HNF3B and TTR. Important for liver genes transcription.



Western blot analysis of extracts from RAW264.7 cells(heat shock treatment), using HNF-6 Ab. The lane on the left was treated with blocking peptide.

IMPORTANT: For western blot, incubate membrane with diluted primary Ab in 5% w/v milk , 1X TBS, 0.1% Tween@20 at 4°C with gentle shaking, overnight.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.